

- / -

1. , /

1.1

	,	
pH		
	100	100~130

(1)

, ,

.

.

(2)

.

,

.

가

.

1

.

(3) pH

가 .

(4)

. 100

가 . ,

100~110 .

, 130

.

가 .

(5)

.

()

.

가

. ,

, 가 .

2

가?

.

.

.

,

가

.

가

1g/l

(80%)

90 , 10~20

.

.

(6)

Levegal PEW

가


,

.

. Oeko tex standard 100

.

1.2



(1) X % :

Y % :

2~4g/l (Levegal PEW, Bayer)

5%

pH 4.5~5 (pH)

(100)

(2) 1~2g/l Taspon WK(, Bayer)

1.3

(1) (100)

Dianix	Yellow AC-E	Resolin	Brilliant Yellow 10GN 400%
	Yellow Brown HRSL-SE 150		Brilliant Yellow 7GL 200%
	Orange UN-SE		Yellow 5GL 200%
	Red FB-E 220		Red FB 200%
	Rubbine UN-SE		Blue FR 02
	Violet 4RS-FS 400		Blue FBL 150%
	Violet 2R-E		Blue Z-GRLN 200%
	Blue GR-E 140		
	Turquoise Blue B-FS 200		
	Navy Blue BG-SE 200		
	Navy Blue BN-E 01		
	Dark Blue B-SE 200		
	Black GS-E 01		

(2) (105~110)

Dianix	Yellow UN-SE 200 New	Resolin	Brilliant Red BLS 200%
	Red UN-SE		Red Violet FBL 200%
	Red H4G-FS		Blue F2GS
	Red HBSL-FS		
	Red 2BSL-FS 150		
	Violet 3R-FS		
	Navy Blue UN-SE200		
	Black BG-FS 200		
	Black RB-FS 200		

(3) ,

Supranol	Yellow 4GL	Isolan	Yellow S-GL
	Orange GSN		Orange S-RL 01
	Red BL 02		Bordeaux S-BLN
	Red 3BW		Navy Blue S-RL
	Red 6BW		Dark Blue 2S-GL 01
	Blue RLW		Brown S-RL
	Blue GLW		Grey S-GL
Telon	Black LDN		
	Black LGN liquid		

2. /

2.1


(Remazol

, DyStar)

1

가

2.2



(1) X% Remazol

Y%

1~3 Avolan REN(, Bayer)

5%

2~4 g/l Levegal PEW

pH 4.5~5(pH)

(100)

(2) 1~2g/l Taspon WK(, Bayer)

2.3

(1) ()

Realan Yellow G gran
Golden Yellow RC gran
Red B gran
Red G gran
Red RC gran
Blue B gran
Blue RC gran
Turquoise G
Navy Blue BG
Black G gran

(2)

(100)

Dianix	Yellow AC-E	Resolin	Brilliant Yellow 10GN 400%
	Yellow Brown HRSL-SE 150		Brilliant Yellow 7GL 200%
	Orange UN-SE		Yellow 5GL 200%
	Red FB-E 220		Red FB 200%
	Rubbing UN-SE		Blue FR 02
	Violet 4RS-FS 400		Blue FBL 150%
	Violet 2R-E		Blue Z-GRLN 200%

Blue GR-E 140

Turquoise Blue B-FS 200

Navy Blue BG-SE 200

Navy Blue BN-E 01

Dark Blue B-SE 200

Black GS-E 01

(105~110)

Dianix	Yellow UN-SE 200 New	Resolin	Brilliant Red BLS 200%
	Red UN-SE		Red Violet FBL 200%
	Red H4G-FS		Blue F2GS
	Red HBSL-FS		
	Red 2BSL-FS 150		
	Violet 3R-FS		
	Navy Blue UN-SE 200		
	Black BG-FS200		
	Black RB-FS 200		